

AEROBIC

BactaServe® Bioculture

Eco-friendly Bacteria ISO 9001:2015 Certified
Proud Recipient-Aqua Foundation Excellence Award 2017



AEROBIC

- ✓ Reduction in commissioning time of ETP & STP
- ✓ Rapid growth of MLSS with high MLVSS Content
- ✓ Acclimatized microbes - Hence higher COD/BOD removal
- ✓ BactaServe Aerobic promotes growth of COD degraders & suppresses growth of Pathogens
- ✓ Shock load stabilization & Performance improvement
- ✓ Reduces BOD & COD and suppresses the foul odour from ETP & STP
- ✓ Reduces Operational costs - No downtime

AREAS OF APPLICATION:

- ✓ ASP - Activated Sludge Process
- ✓ E-ASP - Extended Aeration Process
- ✓ SBR - Sequencing Batch Reactor
- ✓ MBBR - Moving Bed Bio Reactor
- ✓ MBR - Membrane bio Reactor
- ✓ RBC - Rotating Biological Contractor

HOW TO DOSE:

- ✓ Mix 1 kg BactaServe in 20 Litres of feed water. (2Kg in 40 Litres & so on...)
- ✓ Mix it well and dose it as per the detailed dosage schedule
- ✓ Qty required per day depends on wastewater volume & Organic load



Manufactured & Marketed by:
Amalgam Biotech
Division of Amalgam Engineering
Address- Laxman Nagar, Baner, Pune, India

☎ +91-7276020758
✉ sales@amalgambiotech.com
🌐 www.amalgambiotech.com